

\*Tabla de paradas y soporte incluidos en la parte trasera del carro  
\*CPR cardiac board with support included in the back of the trolley

Detalles del producto / Product Details

Superficie superior de trabajo Top work surface	ABS ABS
Empujador Pusher	ABS / ABS
Cajones Drawers	ABS / ABS
Base Base	Compacmel Compacmel
Paredes laterales Side walls	Fibraplast Fibraplast
Perfiles Profiles	Aluminio Anodizado Anodized Aluminium
Guías en U Guides in U	ABS / ABS
Guías Telescópicas Ball bearing slides	Acero Steel material



Accesorios incluidos/ Accessories included

F-NACC-MESA




Mesa extraíble  
Sliding table

F-NACC-02



Empujador de ABS  
Ergonomic ABS pusher

F-NACC-12



Portaseros de 2 ganchos  
IV pole with 2 hooks

F-NACC-18



Soporte para desfibrilador  
Defibrillator support

F-NACC-41



Soporte para bombonas  
Oxygen tank holder

F-NACC-45




Sistema portaetiquetas  
Label holder system

F-NACC-20



Tabla de paradas para  
masaje cardíaco  
CPR cardiac board with  
support

F-ACC-35NOVP



20 uds de cerradura de  
Precinto  
20 pcs of Throw-away seal  
lock